


Joyce Masks Sizing Gauge



INSTRUCTIONS:

- 1: Print out this gauge, ensuring that there is no scaling. Please use the square above to ensure it has printed at 100% scale.
- 2: Cut out the egg shape in the centre.
- 3: Align the top of the sizing gauge with the indent between the top of the nose and the forehead.
- 4: Use the markers to determine the correct size (blue for nasal sizes, green for full face) as shown by the diagrams.

Please note:

This gauge is a reference to aid in sizing and not an exact measurement as other factors may need to be taken into account.

www.EU-PAP.co.uk

intusHEALTHCARE