



A. Measure width of nose:

1"
1.25"
1.5"



Nasal Mask Measurement


- Using the ruler above measure the width of your nose.
- Using your nose measurement, find your mask size in the chart below.

Nose Size	*Advance	Advance Small	Elan®
Small		●	●
Regular	●		●
Large	●		●

*SleepWeaver Advance is one size fits most. If your nose measurement falls under 1.25" try SleepWeaver Advance Small first.

SleepWeaver 3D Neck Sizing Guide


Not actual size
Use cloth tape measure



NOTE: If your measurement falls between sizes on ANY mask, choose the SMALLER mask size.


A. Measure width of nose:

1"
1.25"
1.5"



Full Face Mask Measurements


- Find measurement A. using top ruler.
- Find measurement B. using side ruler.
- Using both measurements, find your mask size in the chart below.



Anew® Size	A. Measurement	B. Measurement
Small	1" to 1.25"	3.25" to 3.5"
Regular	1.25" to 1.5"	3.5" to 3.75"
Large	1.5" & Up	3.75" & Up


B. Measure from bridge of nose to center of chin:

3.25"
3.5"
3.75"




A. Measure width of nose:

25mm
32mm
38mm



51mm



Nasal Mask Measurement


- Using the ruler above measure the width of your nose.
- Using your nose measurement, find your mask size in the chart below.

Nose Size	*Advance	Advance Small	Elan®
Small		●	●
Regular	●		●
Large	●		●

* SleepWeaver Advance is one size fits most. If your nose measurement falls under 1.25" try SleepWeaver Advance Small first.

SleepWeaver 3D Neck Sizing Guide


Not actual size
Use cloth tape measure



NOTE: If your measurement falls between sizes on ANY mask, choose the SMALLER mask size.

A. Measure width of nose:


25mm
32mm
38mm



51mm

Full Face Mask Measurements

- Find measurement A. using top ruler.
- Find measurement B. using side ruler.
- Using both measurements, find your mask size in the chart below.



Anew® Size	A. Measurement	B. Measurement
Small	25 to 32mm	83 to 89mm
Regular	32 to 38mm	89 to 95mm
Large	38mm & Up	95mm & Up

B. Measure from bridge of nose to center of chin:

83mm
89mm
95mm

